




Foley Catheter Holder

- Minimizes catheter movement
- Facilitates catheter traction
- Prevents skin irritation
- Discourages patient tampering
- Latex-free



Features	Benefits
Two inch wide legband, made of exclusive stretch material fits most adults and children without sliding up and down the leg	<ul style="list-style-type: none"> • One size product to fit most patients • Maintains position on the leg, whether ambulatory or bed-ridden • Legband won't narrow or restrict blood flow of deep or superficial veins
Patented Locking Device is Improved Over Earlier Versions	<ul style="list-style-type: none"> • Secures all types of foley catheters, including new generation of silicone and "Teflon" coated • Easier to apply and move foley's position than previous model • Instantly allows clinician to "loop" catheter for greater comfort to prevent meatal movement in the urethra • Allows surgeon to secure foley with traction if desired, (e.g. post-TURP)
Legband Can be Rotated to Move Meatal Pressure Point	<ul style="list-style-type: none"> • Clinician can teach patient to rotate the legband 1" in either direction as needed for better comfort • Nurses can alternate legs on a periodic basis
Improves Skin Care Over Adhesive Tape or Adhesive-Based Holders	<ul style="list-style-type: none"> • Quick, tape-free application • No shaving of male patients required • Eliminates potential for skin breakdown and allergies to tape
Prevents complications that can increase patient length of hospital stay	<ul style="list-style-type: none"> • Reduces nursing labor • Prevents meatal irritation that is a causal factor in nosocomial Urinary Tract Infections (UTI) • Ensures that surgical repair is not compromised by patient tampering • Greater meatal comfort helps prevent patients from pulling out their foley catheter
Latex Free	<ul style="list-style-type: none"> • Prevents patient or care-giver reaction to latex hypersensitivity

Available from



Graykon Scientific
 P.O. Box 200 Mont Albert
 Vic. 3127
 PH: 61 3 98989393
 Fax: 61 3 98989394
 Email: graykon@bigpond.com
 Web: www.users.bigpond.com/graykon

